

HK VAK FILTRA GL

Glass vacuum filtration device

HANSA FLEX

Properties

Connection Hose fitting = 6 mm diameter



Note

Highly viscous media should be diluted with a technical solution (e.g. NHEPTAN-TECH).

Description

Glass vacuum filtration device for visual evaluation of liquids.

Ordering information

Diaphragm filters, vacuum pump and extraction hose must be ordered separately.

Additional information

The oil sample is drawn through the filter and into the filter bottle with the aid of a vacuum pump.
The residue in the filter provides a preliminary indication of the degree of contamination of the oil.
The filtrate bottle can hold up to 1000 ml,
the glass top can hold up to 250 ml,
the filtration area is 12.5 cm² with a filter diameter of 50 mm.

Item

Identification	Weight (kg)
HK VAK-FILTRA-GL	1,7

Accessories

HK VAK MEM Filter for vacuum filtration

Additional elements

HK VAK PUMPE	Vacuum hand pump
PSK	PVC hose, transparent