

# HK VAK FILTRA GL

Glass vacuum filtration device

## Properties

**Connection** Hose fitting = 6 mm diameter



## Note

Highly viscous media should be diluted with a technical solution (e.g. NHEPTAN-TECH).

## Description

Glass vacuum filtration device for visual evaluation of liquids.

## Ordering information

Diaphragm filters, vacuum pump and extraction hose must be ordered separately.

## Additional information

The oil sample is drawn through the filter and into the filter bottle with the aid of a vacuum pump.  
The residue in the filter provides a preliminary indication of the degree of contamination of the oil.  
The filtrate bottle can hold up to 1000 ml,  
the glass top can hold up to 250 ml,  
the filtration area is 12.5 cm<sup>2</sup> with a filter diameter of 50 mm.

## Item

Identification	Weight (kg)
HK VAK-FILTRA-GL	1,7

## Accessories

**HK VAK MEM** Filter for vacuum filtration

## Additional elements

<b>HK VAK PUMPE</b>	Vacuum hand pump
<b>PSK</b>	PVC hose, transparent